


















	MODEL NUMBER	DESCRIPTION	CABLE LENGTH	APPLICATION	PACKAGING	WristOx ₂ [®] Model 3150	PalmSAT [®] 2500/2500A	8500 Series	Model 7500*	Avant [®] 2120	LifeSense [®]	PulseSense [™]
Reusable Soft Sensors 8000S Series												
	8000SS 8000SS-3	Small Soft Sensor 7.5 – 12.5 mm digit 0.3 – 0.5 inches digit	1m 3m	Index, middle or ring finger	Single		✓	✓	✓	✓	✓	✓
	8000SS-WO2	Small Soft Sensor – WristOx₂ 7.5 – 12.5 mm digit 0.3 – 0.5 inches digit	0.3m	Index, middle or ring finger	Single	✓						
	8000SM 8000SM-3	Medium Soft Sensor 10 – 19 mm digit 0.4 – 0.75 inches digit	1m 3m	Index, middle or ring finger	Single		✓	✓	✓	✓	✓	✓
	8000SM-WO2	Medium Soft Sensor – WristOx₂ 10 – 19 mm digit 0.4 – 0.75 inches digit	0.3m	Index, middle or ring finger	Single	✓						
	8000SL 8000SL-3	Large Soft Sensor 12.5 – 25.5 mm digit 0.5 – 1.0 inches digit	1m 3m	Index, middle or ring finger	Single		✓	✓	✓	✓	✓	✓
	8000SL-WO2	Large Soft Sensor – WristOx₂ 12.5 – 25.5 mm digit 0.5 – 1.0 inches digit	0.3m	Index, middle or ring finger	Single	✓						
Reusable Finger Clip, Ear Clip and Reflective Sensors 8000 Series												
	8000AA 8000AA-2 8000AA-3	Finger Clip Sensor – Adult > 30 kg / > 66 lbs	1m 2m 3m	Index, middle or ring finger	Single		✓	✓	✓	✓	✓	✓
	8000AA-WO2	Finger Clip Sensor – WristOx₂ > 30 kg / > 66 lbs	0.3m	Index, middle or ring finger	Single	✓						
	8000AP 8000AP-3	Finger Clip Sensor – Pediatric 8 – 30 kg / 18 – 66 lbs	1m 3m	Index, middle or ring finger	Single		✓	✓	✓	✓	✓	✓
	8000Q2	Ear Clip Sensor > 40 kg / > 88 lbs	1m	Earlobe	Single		✓	✓	✓	✓	✓	✓
	8000R 8000H	Reflectance Sensor > 30 kg / > 66 lbs Sensor Holder 10 caps /20 stickers	1m N/A	Forehead	Single		✓	✓	✓	✓	✓	✓
Reusable Flex Sensors with FlexiWrap 8000J Series												
	8000J 8000J-3 8000JFW	Reusable FlexSensor – Adult > 22 kg / > 44 lbs incl. 25 Adult FlexiWraps	1m 3m N/A	Index, middle or ring finger	Single 24 per pack		✓	✓	✓	✓	✓	✓
	8000J-WO2 8000JFW	Reusable FlexSensor – WristOx₂ > 30 kg / > 66 lbs incl. 25 Adult FlexiWraps	0.3m N/A	Index, middle or ring finger	Single 24 per pack	✓						
	8008J 8008JFW	Reusable FlexSensor – Infant 2 – 22 kg / 4.4 – 44 lbs incl. 25 Infant FlexiWraps	1m N/A	Toe	Single 24 per pack		✓	✓	✓	✓	✓	✓
	8001J 8001JFW	Reusable FlexSensor – Neonatal > 2 kg / 4.4 lbs incl. 25 Neonatal FlexiWraps	1m N/A	Side of foot	Single 24 per pack		✓	✓	✓	✓	✓	✓
Single-Patient-Use Cloth Sensors 6000C Series												
	6000CA	Cloth Sensor – Adult > 30 kg / > 66 lbs	1m	Index, middle or ring finger	Single 24 per pack		✓	✓	✓	✓	✓	✓
	6000CP	Cloth Sensor – Pediatric 10 – 40 kg / 22 – 88 lbs	1m	Index, middle or ring finger	Single 24 per pack		✓	✓	✓	✓	✓	✓
	6000CI	Cloth Sensor – Infant 2 – 20 kg / 4.4 – 44 lbs	1m	Big toe	Single 24 per pack		✓	✓	✓	✓	✓	✓
	6000CN	Cloth Sensor – Neonatal/Adult Neonatal: < 2 kg / < 4.4 lbs Adult: > 30 kg / > 66 lbs	1m	Foot	Single 24 per pack		✓	✓	✓	✓	✓	✓
Single-Patient-Use Flexi-Form[®] III Sensors 7000 Series												
	7000A	Flexi-Form Sensor – Adult > 30 kg / > 66 lbs	1m	Index, middle or ring finger	Single 24 per pack		✓	✓	✓	✓	✓	✓
	7000P	Flexi-Form Sensor – Pediatric 10 – 40 kg / 22 – 88 lbs	1m	Index, middle or ring finger	Single 24 per pack		✓	✓	✓	✓	✓	✓
	7000I	Flexi-Form Sensor – Infant 2 – 20 kg / 4.4 – 44 lbs	1m	Big toe	Single 24 per pack		✓	✓	✓	✓	✓	✓
	7000N	Flexi-Form Sensor – Neonatal/Adult Neonatal: < 2 kg / < 4.4 lbs Adult: > 30 kg / > 66 lbs	1m	Foot	Single 24 per pack		✓	✓	✓	✓	✓	✓
Single-Patient-Use Durafoam Sensors 6500 Series												
	6500SA	Durafoam Sensor – Standard 5 – 21 mm/0.19 – 0.82 inches digit > 30 kg / > 66 lbs	1m	Finger	5 per pack		✓	✓	✓	✓	✓	✓
	6500MA	Durafoam Sensor – Small 5 – 21 mm/0.19 – 0.82 inches digit > 30 kg / > 66 lbs	1m	Finger	5 per pack		✓	✓	✓	✓	✓	✓
Accessories												
	UNI-EXT-1 UNI-EXT-3 UNI-EXT-6 UNI-EXT-9	Extension Cable	1m 3m 6m 9m		Single		✓	✓	✓	✓	✓	✓
	UNI-RA-3 UNI-RA-0	Right Angle Connector	3m 15cm		Single		✓	✓	✓	✓	✓	✓
	31501	Adapter cable to 9-pin (D-Connector). Connects PureLight [®] sensors to the WristOx ₂ Model 3150	12cm		Single	✓						
	3150WI	Adapter cable for Avant 4100 pulse oximeter and WristOx 3100 sensor	12cm		Single							

*Also available is the Model 7500FO for use in MRI environments utilizing the 8000FI and 8000FC cable and sensor.