







TALKING + UPPER ARM BLOOD PRESSURE MONITOR



KEY FEATURES:

-  TriCheck
-  Coaching
-  Irregular Heartbeat Detection
-  Blood Pressure Rating Indicator
-  90 Reading Memory
-  5 Year Warranty



CLINICALLY
VALIDATED

A premium talking blood pressure monitor with advanced inflation technology that provides audible readings and instructions in multiple languages.

Trusted Accuracy

- Clinically validated for accuracy
- TriCheck™ mode averages 3 consecutive readings per medically recommended guidelines
- Irregular Heartbeat (IHB) alerts the user to irregular heartbeat patterns
- Pressure Rating Indicator classifies blood pressure in accordance with the medically approved standards
- Cuff fit error prevents inaccurate measurements by alerting user when the cuff is too loose
- Body movement sensor guides user to remain still during measurement

User Friendly

- Talking multi-lingual device gives verbal instructions and results
- Easy, one button operation for quick and accurate measurement
- Stores 90 measurements in memory

- Easy to use medium cuff (9.0-14.6") that conforms to the user's arm to provide a comfortable experience
- Adjustable pressure settings
- Advanced inflation and deflation technology providing faster and more comfortable readings
- Clock display with date and time stamp for logging and trending measurements
- Large display is easy to read
- Includes an AC adapter
- Carrying case included

The A&D Promise

- Professional grade performance you can trust
- 5 year warranty on monitor, 2 year warranty on cuff
- Cuffs are not made with natural rubber latex
- Battery life indicator
- Download the A&D Heart Track App for Free to manually log and electronically track measurements (available on iOS and Android)

Technical Specification

Measurement Method	Oscillometric
Memory	90 measurements
Cuff Size	9.0-14.6" (23-37 cm)
Measurement Range	Pressure: 0-280 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-200 beats/min
Measurement Accuracy	Pressure: ±3 mmHg Pulse: ±5 %
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-233)
Battery Life	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)
Monitor Dimensions	5.5 x 4.1 x 2.3" (14 x 10 x 6 cm)
Monitor Weight	10.6 oz (300g) without batteries
Display Size	2.4 x 2.8" (6.1 x 7.1 cm)
Clinical Test	ANSI / AAMI SP-10 1992
Operating Conditions	+10°C to +40°C (50°F to 104°F) / 15% to 85% RH / 800 to 1060 hPa
Storage/Transport Conditions	-20°C to +60°C (-4°F to 140°F) / 10% to 95% RH / 700 to 1060 hPa

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
UA-1030T Individual Unit (Box)	7.1 x 4.0 x 8.5" (18 x 10 x 22 cm)	2.1 lbs (953 g)	0.14 ft³ (0.004 m³)	0-93764-60293-1
UA-1030T Case (10 units)	21 x 15 x 9.5" (53 x 38 x 24 cm)	22.6 lbs (10.2 kg)	1.73 ft³ (0.05 m³)	1-00-93764-60293-8

What is In The Box

- 1 Blood Pressure Monitor
- 1 Blood Pressure Cuff (9.0-14.6")
- 1 Manual
- 1 Quick Start Guide
- 1 Carrying Case
- 4 AA Batteries
- AC Adapter